

Project1

February 21, 2016

```
In [4]: from QTClock import QTClock
        from ECGFigure import ECGFigure
        from HRDerivClock import HRDerivClock
        from ECGClock import ECGClock
        import matplotlib.pyplot as plt
        import matplotlib.image as mpimg
        %matplotlib inline
        import pandas as pd
        from scipy import signal
        import numpy as np
        import math

In [5]: def smoothhr(inputdata,outputdata,cut = 10,shifttime=0,filtering = 10):
        th =pd.read_csv(inputdata, sep=',')

        #leads = len(th['lead'].unique())
        th = th[th.hearttrate >= 40] #throw hearttrate smaller than 40 or bigger than 160
        th = th[th.hearttrate <= 160]

        #th['HeartRate1']= pd.rolling_median(th['HeartRate'], window=9, center=True).fillna(method=
        total = len(th['hearttrate'])

        time = th.Time.tolist()

        h = []
        t = []
        for i in range(0,total-filtering,filtering):
            strt = i
            end = (i+filtering)

            ti= time[strt]
            hi = th['hearttrate'][strt:end].median()

            t.append(ti)
            h.append(hi)

        t = pd.to_datetime(t, dayfirst=True)
        t = t + pd.DateOffset(hours=shifttime)

        result = pd.DataFrame(h,t)

        result = result.ix[cut:]
```

```

result = result.ix[:~cut]

result.to_csv(outputdata,header = False) #remove the first one

def cuberoot(x):
    return math.pow(x,1/3)

def smoothqtcf(inputdata,outputdata,cut = 10,shifttime=0,filtering = 10):
    th =pd.read_csv(inputdata, sep=',')

    th = th[th.heartrate >= 40] #throw heartrate smaller than 40 or bigger than 160
    th = th[th.heartrate <= 160]
    th = th[th.qti >= 210] #throw sqti smaller than 700 or bigger than 210
    th = th[th.qti <= 700]

    total = len(th['qti'])

    th['rri'] = th['rri'].astype('float64')
    th['qti'] = th['qti'].astype('float64')

    #th['rri_s'] = th['rri'].apply(np.sqrt)
    th['rri_s3'] = th['rri'].apply(cuberoot)
    th['qtcF'] = th['qti']/th['rri_s3']

    time = th.Time.tolist()

    f = []
    t = []
    for i in range(0,total-filtering,filtering):
        strt = i
        end = (i+filtering)
        ti = time[strt]
        fi = th['qtcF'][strt:end].median()

        t.append(ti)
        f.append(fi)

    t = pd.to_datetime(t, dayfirst=True)
    t = t + pd.DateOffset(hours=shifttime)

    #f = th.qtcF.tolist()
    resultf = pd.DataFrame(f,t)

    resultf = resultf.ix[cut:]
    resultf = resultf.ix[:~cut]

    resultf.to_csv(outputdata,header = False)

In [3]: smoothqtcf('RRQT/9004_rrqt.csv','qtcf/9004_qtcf.csv')

In [4]: my_clock = QTClock('QTcf patient9004 health male 28yr')

my_clock.add_percentile_range('./normal_ranges/QTcF_healthy_male.csv',
                             lower=7, upper=93, color='g', alpha=0.2, label='healthy')

```

```


my_clock.add_recording('qtcf/9004_qtcf.csv', label='9004.ecg', filtering = 10)

# >500ms will be highlighted red:
my_clock.add_danger_range(500)

my_clock.add_legend()

my_clock.save('photos/9004_qtcf.png')

```



```

In [5]: #figure 1
smoothqtcf('RRQT/206_rrqt.csv', 'qtcf/206_qtcf.csv', cut = 10, shifttime=15, filtering = 10)

```

```

In [6]: #figure 1

```

```

my_clock = QTClock('QTcf patient206 1-year-old female LQT2')

my_clock.add_percentile_range('./normal_ranges/QTcF_healthy_female.csv',
                              lower=7, upper=93, color='g', alpha=0.2, label='healthy')

#my_clock.add_percentile_range('./normal_ranges/QTcF_LQT1_female.csv', color='b', alpha=0.4, la

```

```

my_clock.add_percentile_range('./normal_ranges/QTcf_LQT2_female.csv', color='y', alpha=0.4, lab


my_clock.add_recording('qtcf/206_qtcf.csv', label='206.ecg', filtering = 10)

# >500ms will be highlighted red:
my_clock.add_danger_range(500)

my_clock.add_legend()

my_clock.save('photos/206_qtcf_inpaper.png')

```



```

In [7]: #figure 2
smoothqtcf('RRQT/288_rrqt.csv', 'qtcf/qtcf_288.csv', shifttime=17)

```

```

In [8]: #figure 2
#smoothqtcf('RRQT/288_rrqt.csv', 'qtcf/qtcf_288.csv', shifttime=17)

```

```

my_clock = QTClock('QTcf (ms) for 32-year-old female LQT1', color_cycle = ['b', 'c'])

```

```

# Add two recordings to the plot. Filtering is disabled because this data was
# already filtered.

# Show two 'tiers' of healthy QTc ranges: IQR in darker green, and a wider
# percentile range in lighter green around it. IQR is darker because the
# regions overlap. We only label one of them, because we don't want redundant
# entries in the legend.
my_clock.add_percentile_range('./normal_ranges/QTcF_healthy_female.csv',
                             lower=16, upper=84, color='g', alpha=0.2, label='healthyFemale')

my_clock.add_percentile_range('./normal_ranges/QTcF_LQT1_female.csv',
                             lower=16, upper=84, color='y', alpha=0.2, label='LQT1')


my_clock.add_recording('qtcf/qtcf_288.csv', label='unfiltered')
my_clock.add_recording('qtcf/qtcf_288.csv', label='filtered', filtering = 10)

# >500ms will be highlighted red:
my_clock.add_danger_range(500)

my_clock.add_legend()

my_clock.save('qtcf/qtcf_288_inpaper.png')

```



```
In [9]: smoothqtcf('RRQT/109_rrqt.csv','qtcf/109_qtcf.csv',shifttime=10,filtering = 10)
smoothqtcf('RRQT/10_rrqt.csv','qtcf/10_qtcf.csv',shifttime=10,filtering = 10)
smoothqtcf('RRQT/8_rrqt.csv','qtcf/8_qtcf.csv',shifttime=10,filtering = 10)
smoothqtcf('RRQT/7_rrqt.csv','qtcf/7_qtcf.csv',shifttime=10,filtering = 10)
smoothqtcf('RRQT/9_rrqt.csv','qtcf/9_qtcf.csv',shifttime=10,filtering = 10)
```

```
In [10]: #figure 4
```

```
my_fig = ECGFigure(title='QTcF (ms) for THEW ID 111',ncols=2)
```

```
left_clock = QTClock('Consistent QTc', parent_figure = my_fig,subplot=1)
```

```
right_clock= QTClock('Inconsistent QTc',parent_figure = my_fig, subplot =2)
```

```
left_clock.add_percentile_range('./normal_ranges/QTcF_healthy_male.csv',color='g', alpha=0.4, l
```

```
right_clock.add_percentile_range('./normal_ranges/QTcF_healthy_male.csv', color='g', alpha=0.4
```

```
left_clock.add_percentile_range('./normal_ranges/QTcF_LQT2_male.csv', color='y', alpha=0.4, l
```

```
right_clock.add_percentile_range('./normal_ranges/QTcF_LQT2_male.csv', color='y', alpha=0.4, l
```

```


left_clock.add_recording('qtcf/109_qtcf.csv', label='age 3', filtering = 20, color = 'b')
left_clock.add_recording('qtcf/10_qtcf.csv', label='age 4', filtering = 20, color = 'm')
left_clock.add_recording('qtcf/8_qtcf.csv', label='age 6', filtering = 20, color = 'c')

right_clock.add_recording('qtcf/7_qtcf.csv', label='age 5', filtering=20, color='k')
right_clock.add_recording('qtcf/8_qtcf.csv', label='age 6', filtering=20, color='c')
right_clock.add_recording('qtcf/9_qtcf.csv', label='age 7', filtering=20, color='y')

left_clock.add_legend()
right_clock.add_legend()

my_fig.save('photos/qtcf_figure4_inpaper.png')

```



```

In [11]: smoothqtcf('RRQT/268_rrqt.csv', 'qtcf/268_qtcf.csv', filtering = 10)
In [12]: smoothqtcf('RRQT/377_rrqt.csv', 'qtcf/377_qtcf.csv', filtering = 10)
In [13]: smoothqtcf('RRQT/427_rrqt.csv', 'qtcf/427_qtcf.csv', filtering = 10)
In [14]: my_clock = QTClock('QTcf LQT1 female 14yrs vs 28yrs vs 3month', color_cycle = ['b', 'c', 'm'])

my_clock.add_percentile_range('./normal_ranges/QTcF_healthy_female.csv',
                              lower=16, upper=84, color='g', alpha=0.2, label='healthyFemale')

my_clock.add_percentile_range('./normal_ranges/QTcF_LQT1_female.csv', color='y', alpha=0.2, label='LQT1Female')


my_clock.add_recording('qtcf/268_qtcf.csv', label='14yrs patient 268', filtering = 10)
my_clock.add_recording('qtcf/377_qtcf.csv', label='28yrs patient 377', filtering = 10)
my_clock.add_recording('qtcf/427_qtcf.csv', label='3mon patient 427', filtering = 10)

```

```
# >500ms will be highlighted red:
my_clock.add_danger_range(500)

my_clock.add_legend()

my_clock.save('photos/268vs377vs427.png')
```



```
In [15]: smoothqtcf('RRQT/208_rrqt.csv','qtcf/208_qtcf.csv',filtering = 20)
In [16]: smoothqtcf('RRQT/183_rrqt.csv','qtcf/183_qtcf.csv',filtering = 20)
In [17]: smoothqtcf('RRQT/170_rrqt.csv','qtcf/170_qtcf.csv',filtering = 20)
In [18]: my_clock = QTClock('QTcf LQT2 male 38yrs vs 10yrs vs 72 yrs',color_cycle = ['b','c','r'])

my_clock.add_percentile_range('./normal_ranges/QTcF_healthy_male.csv',
                              lower=16, upper=84, color='g', alpha=0.2, label='healthyFemale')

my_clock.add_percentile_range('./normal_ranges/QTcF_LQT2_male.csv',
```



```
lower=16, upper=84, color='y', alpha=0.2, label='LQT2')
```


```
my_clock.add_recording('qtcf/208_qtcf.csv', label='38yrs patient 208', filtering = 10)
my_clock.add_recording('qtcf/183_qtcf.csv', label='10yrs patient 183', filtering = 10)
my_clock.add_recording('qtcf/170_qtcf.csv', label='72yrs patient 170', filtering = 10)
```

```
# >500ms will be highlighted red:
```

```
my_clock.add_danger_range(500)
```

```
my_clock.add_legend()
```

```
my_clock.save('photos/208vs183vs170.png')
```



```
In [19]: smoothqtcf('RRQT/273_rrqt.csv', 'qtcf/273_qtcf.csv', filtering = 20)
```

```
In [20]: my_clock = QTClock('QTcf LQT1 6yr male vs LQT2 6yr male ', color_cycle = ['b', 'c'])
```

```
my_clock.add_percentile_range('./normal_ranges/QTcF_healthy_male.csv',
```

```

lower=16, upper=84, color='g', alpha=0.2, label='healthyFemale')

my_clock.add_percentile_range('./normal_ranges/QTcf_LQT1_male.csv',
                              lower=16, upper=84, color='m', alpha=0.2, label='LQT1')

my_clock.add_percentile_range('./normal_ranges/QTcf_LQT2_male.csv',
                              lower=16, upper=84, color='y', alpha=0.2, label='LQT2')


my_clock.add_recording('qtcf/8_qtcf.csv', label='LQT2 6yr male patient 8', filtering = 10)
my_clock.add_recording('qtcf/273_qtcf.csv', label='LQT1 6yr male patient 273', filtering = 10)

# >500ms will be highlighted red:
my_clock.add_danger_range(500)

my_clock.add_legend()

my_clock.save('photos/8vs273.png')

```



```

In [21]: smoothhr('RRQT/9004_rrqt.csv', 'hr/9004_hr.csv', filtering = 10, shifttime = 10)
In [22]: #change the title below
my_fig = ECGFigure(nrows=1, ncols=2, title='health heart rate 9004.ecg')

```


```

hr_clock = ECGClock('HR(bpm)', parent_figure=my_fig, subplot=1,min_rad = 40,max_rad = 105)
hr_clock.add_recording('hr/9004_hr.csv',filtering = 10,label='9004')

#change the range for different type
hr_clock.add_percentile_range('./normal_ranges/HR_healthy_male.csv',color = 'g',
                             alpha = 0.4, label = "healthM")

dhr_clock = HRDerivClock('dHR/dt (percent/second)', parent_figure = my_fig, subplot =2,
                          min_rad = -1, max_rad = 1, tick_spacing = 1.0)
dhr_clock.add_recording('hr/9004_hr.csv',percent = True)
hr_clock.add_legend()
my_fig.save('photos/9004hr_inpaper.png')

```



```

In [23]: smoothhr('RRQT/206_rrqt.csv','hr/206_hr.csv',filtering = 20,shifttime = 0)

```

```

In [24]: #change the title below

```

```

my_fig = ECGFigure(nrows=1, ncols=2, title='hr LQT2 1yr female')


hr_clock = ECGClock('HR(bpm)', parent_figure=my_fig, subplot=1,min_rad = 40,max_rad = 160)
hr_clock.add_recording('hr/206_hr.csv',filtering = 10,label='patient 206')

#change the range for different type

hr_clock.add_percentile_range('./normal_ranges/HR_LQT2_female.csv',color = 'r',
                             alpha = 0.4, label = "LQT2Female")

dhr_clock = HRDerivClock('dHR/dt (percent/second)', parent_figure = my_fig, subplot =2,
                          min_rad = -1, max_rad = 1, tick_spacing = 1.0)
dhr_clock.add_recording('hr/206_hr.csv',percent = True)
hr_clock.add_legend()
my_fig.save('photos/206hr.png')

```



```
In [25]: smoothhr('RRQT/288_rrqt.csv','hr/288_hr.csv',filtering = 20,shifttime = 0)
```

```
In [26]: smoothhr('RRQT/268_rrqt.csv','hr/268_hr.csv',filtering = 20,shifttime = 0)
```

```
In [27]: smoothhr('RRQT/377_rrqt.csv','hr/377_hr.csv',filtering = 20,shifttime = 0)
```

```
In [28]: smoothhr('RRQT/427_rrqt.csv','hr/427_hr.csv',filtering = 20,shifttime = 0)
```

```
In [7]: hr_clock = ECGClock('LQT1 female 32yrs vs 14yrs vs 28 yrs vs 3mon', min_rad = 30,max_rad = 170)
hr_clock.add_recording('hr/288_hr.csv',filtering = 10,label='32yrs patient 288',color = 'b')
hr_clock.add_recording('hr/268_hr.csv',filtering = 10,label='14yrs patient 268',color = 'c')
```


```
hr_clock.add_recording('hr/377_hr.csv',filtering = 10,label='28yrs patient 377',color = 'y')
```

```
hr_clock.add_recording('hr/427_hr.csv',filtering = 10,label='3mon patient 427',color = 'm')
```

```
#change the range for different type
```

```
hr_clock.add_legend()
```

```
hr_clock.save('photos/hr_LQT1_female.png')
```




```
In [30]: smoothhr('RRQT/109_rrqt.csv','hr/109_hr.csv',filtering = 20,shifttime = 11)
In [31]: smoothhr('RRQT/10_rrqt.csv','hr/10_hr.csv',filtering = 20,shifttime = 11)
In [32]: smoothhr('RRQT/7_rrqt.csv','hr/7_hr.csv',filtering = 20,shifttime = 11)
In [33]: smoothhr('RRQT/8_rrqt.csv','hr/8_hr.csv',filtering = 20,shifttime = 11)
In [34]: smoothhr('RRQT/9_rrqt.csv','hr/9_hr.csv',filtering = 20,shifttime = 11)
In [35]: hr_clock = ECGClock('HR same patient LQT2 male', min_rad = 40,max_rad = 140)
hr_clock.add_recording('hr/109_hr.csv',filtering = 10,label='3yr patientid 109',color = 'b')
hr_clock.add_recording('hr/10_hr.csv',filtering = 10,label='5yr patientidn 7',color = 'c')

hr_clock.add_recording('hr/7_hr.csv',filtering = 10,label='7yr patientid 9',color = 'y')
hr_clock.add_recording('hr/8_hr.csv',filtering = 10,label='6yr patientid 8',color = 'm')
hr_clock.add_recording('hr/10_hr.csv',filtering = 10,label='4yr patientid 10',color = 'g')

hr_clock.add_percentile_range('./normal_ranges/HR_LQT2_male.csv',color = 'r',
                             alpha = 0.4, label = "LQT2 male")
#change the range for different type
hr_clock.add_legend()
```

```
hr_clock.save('photos/hr_LQT2_male_samepatient.png')
```




```
In [36]: smoothhr('RRQT/273_rrqt.csv','hr/273_hr.csv',filtering = 20,shifttime = 12)
```

```
In [9]: hr_clock = ECGClock('6yrs male LQT2 vs LQT1', min_rad = 40,max_rad = 130)
hr_clock.add_recording('hr/8_hr.csv',filtering = 10,label='LQT2 patientid 8',color = 'b')
hr_clock.add_recording('hr/273_hr.csv',filtering = 10,label='LQT1 patientidn 273',color = 'c')

hr_clock.add_percentile_range('./normal_ranges/HR_LQT1_male.csv',color = 'y',
                              alpha = 0.4, label = "LQT1 male range")
hr_clock.add_percentile_range('./normal_ranges/HR_LQT2_male.csv',color = 'r',
                              alpha = 0.4, label = "LQT2 male range")
#change the range for different type
hr_clock.add_legend()

hr_clock.save('photos/hr_6yrs_male_LQT2_vs_LQT1.png')
```



```
In [38]: smoothhr('RRQT/208_rrqt.csv','hr/208_hr.csv',filtering = 20,shifttime = 12)


In [39]: smoothhr('RRQT/183_rrqt.csv','hr/183_hr.csv',filtering = 20,shifttime = 12)

In [40]: smoothhr('RRQT/170_rrqt.csv','hr/170_hr.csv',filtering = 20,shifttime = 12)

In [8]: hr_clock = ECGClock('HR LQT2 male 38yrs vs 10yrs vs 72yrs', min_rad = 40,max_rad = 140)
hr_clock.add_recording('hr/208_hr.csv',filtering = 10,label='38yr patientid 208',color = 'b')
hr_clock.add_recording('hr/183_hr.csv',filtering = 10,label='10yr patientidn 183',color = 'c')
hr_clock.add_recording('hr/170_hr.csv',filtering = 10,label='72yr patientid 170',color = 'y')

hr_clock.add_percentile_range('./normal_ranges/HR_LQT2_male.csv',color = 'r',
                             alpha = 0.4, label = "LQT2 male")
#change the range for different type
hr_clock.add_legend()

hr_clock.save('photos/hr_LQT2_male_differences.png')
```



In []:

In []: